



The Endomed 482 is the electrotherapy device for the demanding therapist.

Electrotherapy has many possibilities and applications. The deployment of pain control and muscle stimulation has a long history and is well substantiated scientifically. The use of electrotherapy is not invasive, easy to apply and will not burden the patient. Electrotherapy is, therefore, often an excellent alternative to medication or other interventions.

The Endomed 482 is characterised by its convenient use because of its touch screen interface and the extensive protocol library (no fewer than 42 items and evidence based) as is the case with the other devices from the Enraf-Nonius 4 series.

The Endomed 482 can be expanded with a Vacotron 460. This integrated unit ensures that applying electrotherapy by using vacuum electrodes becomes very easy.

The Endomed 482: the complete, easy-to-use, fast and safe device for electrotherapy!



ELECTROTHERAPY

Equipment for pain management, muscle stimulation and/or diagnostics. Electrotherapy applications for fast pain relief without side effects. Muscles can be stimulated to improve the function of soft tissue. Electrotherapy can be used to measure differentiated tissue response to stimulation, to provide a diagnosis about tissue quality.

0498732_40 EN 01052017

1498920 ENDOMED 482



TECHNICAL SPECIFICATIONS

Power supply 100-240 V \pm 10% / 50 / 60 Hz Dimensions (WxDxH) Device 24 x 32 x 12 cm Device with base 24 x 30.5 x 18.2 cm Device with Vacotron 24 x 30.5 x 21.6 cm Weight 2 kg Weight, including the optional battery 3 kg

ELECTROTHERAPY

Current channels Current types Pre-programmed treatment suggestions 42 - evidence based 1000+ Programmable positions

2 (100% independent) 26 (per channel)

OTHER PRODUCTS THAT MIGHT INTEREST YOU

VACOTRON 460

The Sonopuls 492 can be extended with a vacuum module, the Vacotron 460. A vacuum unit for easy application of vacuum electrodes for combination therapy and electrotherapy.

EN-CAR U3

The EN-Car U3 is your companion for any piece of Enraf-Nonius equipment. For laser, ultrasound, biofeedback, shockwave or electrotherapy, with the EN-Car U3 you have it all close at hand.

ORDERING DATA

1498950 Vacotron 460 1468960 EN-Car U3

Also see the seperate product leaflets at WWW.ENRAF-NONIUS.COM/ELECTROTHERAPY



STANDARD ACCESSORIES ENDOMED 482

	1498757	4-Series Operating Instructions (cd rom)
	1498756	4-Series Information Booklet
	3444290	Mains cable 230V-EUR
	1498010	Device base (inclination support) (not for Vacotron 460)
	3444021	Strap 250x3 cm
	3444020	Strap 100x3 cm
2x	3444129	Rubber electrodes 6x8 cm, 2 mm female, set of 2
	1460266	Moist pads for rubber electrodes 6x8 cm, set of 4
2x	3444211	Patient cable 2-core & 2 mm male plugs - black,
		with colour clips







ENDOMED 482

• • • • • •



42 Evidence-Based protocols for electrotherapy. Complete including the display of the anatomic localisations and the treatment instructions.





REMOTE CONTROL (OPTIONAL)

The patient can be involved in his or her own treatment by using the optional "remote control". Especially when using electrotherapy, the patient will have the possibility to control the therapy.

BATTERY-OPERATED (OPTIONAL) AND, THERE-FORE, MOBILE

In addition to being able to use the 4-series unit connected to a mains power supply, the device can also be used with a battery-operated power supply. This will ensure that you can go everywhere with your unit, regardless of whether this is at the patient's home or at his or her hospital bed.








MADE BY ENRAF-NONIUS

Made by Enraf-Nonius and, therefore, safe, reliable and durable. A device made by Enraf-Nonius, THE supplier for the physiotherapist worldwide. Where quality, reliability and durability come as standard.



CRYSTAL CLEAR TFT COLOUR SCREEN





Easy to use by means of a clear menu via a touch screen interface.

TOUCH SCREEN USER INTERFACE

Extremely intuitive and as easy as a state-of-the-art smart phone.

COMPACT AND STYLISH

An extremely compact device, small and light-weight for easy transportation. At the same time stylish and beautiful, a device that is worth-while showing!





PROGRAMMED FOR SPEED

ENDOMED 482

• • • • • •

The menu structure is optimized to speed. In most cases you can start a pathology-specific treatment in just 3 steps! Easy, quick and simple.




ENRAF NONIUS



PROGRAMMABLE

In addition to the standard treatment protocols, you can also save your personal protocols (including sequential applications) based on your own guidelines.



ACCESSORIES ENDOMED 482



	DES AND MOIST PADS				
ADHESIVI	ELECTRODES				
3444222	Adhesive electrodes Ø 2,2 cm, 2 mm female, 10 sheets of 8 (also for EMG)	3444222	3444056	3444135	3444057
3444056	EN-Trode Ø 3,2 cm, 2 mm female,10 sheets of 4	5111222 1	5111050	541155	544057
3444135	EN-Trode Ø 5,0 cm, 2 mm female, 10 sheets of 4				
3444057	EN-Trode 5x5 cm, 2 mm female, 10 sheets of 4				
3444058	EN-Trode 5x9 cm, 2 mm female, 10 sheets of 4				
RUBBER ELECTRODES		3444058			
3444128	Rubber electrodes 4x6 cm, 2 mm female, set of 2.				
3444129	Rubber electrodes 6x8 cm, 2 mm female, set of 2			-	
3444130	Rubber electrodes 8x12 cm, 2 mm female, set of 2			On	
MOIST PA	DS FOR RUBBER ELECTRODES		_		
1460273	Moist pads for rubber electrodes 4x6 cm, set of 4	3444128	3444129	3444130	
1460266	Moist pads for rubber electrodes 6x8 cm, set of 4				
1460275	Moist pads for rubber electrodes 8x12 cm, set of 4				
FIXATION	STRAPS				
3444020	Strap 100x3 cm		^	1	
3444021	Strap 250x3 cm	1460273	1460266	1460275	
3444022	Strap 100x5 cm				
3444023	Strap 250x5 cm				
POINT ELI	CTRODES				
1480800	Electrode holder for 1480801	Τ	Т	т 3444022	Т 3444023
1480801	Point electrode 5 mm, 4 mm female (in combination with adapters 2523523 / 2523524)	3444020	3444021	5444022	544025
		1480800	1480801		
ADAPTER	S				
2523524	Adapter plug, 2 mm female, 4 mm male, red	100			
2523523	Adapter plug, 2 mm female, 4 mm male, black	2523524	2523523		
PATIENT	CABLES	~ -			
3444211	Patient cable 2-core & 2 mm male plug - black, with coloured clips	3444211			
REMOTE	CONTROLS				
1498800	Remote control, for 4-series				
1150000		1498800			
CASES AN	D BAGS	1 And			
3444675	Carrierbag 4-series	3444675			
AC ADAP	TERS				
AC ADAP					

2501016

• • • • • •

ENRAF-NONIUS